

04-2015 / 9308960/002/CB/B / This document is not legally binding. bioMérieux reserves the right to modify specifications without notice. / bioMérieux S.A. the blue logo, API 20 E, API 20 A, APIMEB and API are used, pending and/or registered trademarks belonging to bioMérieux S.A. or one of its subsidiaries. / Any other name or trademark is the property of its respective owner/ bioMérieux SA 673 620 399 RCS Lyon / Printed in France / RCS Lyon B 398 160 242

bioMérieux S.A.
69280 Marcy l'Etoile
France
Tel. : 33 (0)4 78 87 20 00
Fax : 33 (0)4 78 87 20 90
www.biomerieux.com



API & ID 32 IDENTIFICATION DATABASES |



INTRODUCTION

The API[®], ID 32 and rapid ID 32 database update takes into account:

- the evolution of international taxonomy
- the description of the new bacterial species,
- newly acquired bacteriology data (new profiles for bacterial strains which have an impact on performance data)

As a result of the update, the APIWEB™ software version has changed from version 1.2.1 to version 1.3.0

The API and ID 32 databases have again been updated

Twenty-two of the twenty-three identification databases have been revised, taking account the biochemical profiles of over 56,277 strains. Today, 697 species of bacteria and yeasts can be identified, including 14 new species and 50 that have been assigned new names.

WHAT'S CHANGED IN THE DATABASES?

The changes made can be broken down as follows:

A number of new species have been added to the database (including both entirely new species and others added on the basis of new results).

Certain bacterial species have been deleted due to more stringent criteria. Certain rare species which are not sufficiently studied have been removed from the database.

The names of certain species have been changed to follow modifications in the bacterial taxonomy as officially described in the International Journal of Systematic and Evolutionary Microbiology.

Notes have been revised to reflect the changes in names and the species added and deleted.


Percentages and performances have been altered to reflect variations observed in the profiles analyzed as the database was revised.

Additional tests were modified to reflect the new reference information available.

Database	Version number		Changes to thesaurus		Changes to database	
	old	new	Taxons	Notes	Identification	Additional tests
API® 20 E	v 4.1	v 5.0	X	X	X	X
RapiD 20 E™	v 3.1	v 3.2	X	-	-	-
API® 10 S	v 3.1	v 4.0	X	X	X	-
API® 20 NE	v 7.0	v 8.0	X	X	X	-
API® STAPH	v 4.1	v 5.0	-	-	X	-
API® 20 STREP	v 7.0	v 8.0	X	X	X	X
API® 20 C AUX	v 4.0	v 5.0	X	X	X	-
API® CANDIDA	v 2.1	v 2.2	X	X	-	-
API® 20 A	v 4.0	v 5.0	X	X	X	X
API® CORYNE	v 3.0	v 4.0	X	X	X	-
API® CAMPY	v 2.1	v 3.0	X	-	X	-
API® LISTERIA	v 1.2	v 2.0	-	X	X	-
API® NH	v 3.0	v 4.0	X	X	X	-
API® 50 CHB	v 4.0	v 4.1	X	-	-	-
API® 50 CHE	v 3.1	v 3.2	X	-	-	-
API® 50 CHL	v 5.1	v 5.2	X	-	-	-
ID 32 E	v 3.0	v 4.0	X	-	X	X
rapid ID 32 E	v 3.1	v 4.0	X	-	X	X
ID 32 STAPH	v 2.1	v 3.0	X	X	X	X
rapid ID 32 STREP	v 3.0	v 4.0	X	X	X	X
ID 32 C	v 3.0	v 4.0	X	X	X	X
rapid ID 32 A	v 3.2	v 3.3	X	-	-	-


ORIENTATION TESTS

GRAM⁻ BACILLI




ORIENTATION TESTS

OTHER GRAM⁻ BACILLI




ORIENTATION TESTS

GRAM⁺ BACILLI



ORIENTATION TESTS

COCCI



ORIENTATION TESTS

YEASTS



Assimilation tests

api 20C AUX
ID 32 C

Candida
Cryptococcus
Geotrichum
Kloeckera
Pichia
Rhodoturela
Trichosporon

Sugar acidification or enzymatic tests

api CANDIDA

Candida
Cryptococcus
Saccharomyces
Trichosporon

SPECIES IDENTIFIABLE

BY THE VARIOUS IDENTIFICATION SYSTEMS

API® 20 E	Gram-negative bacilli
API® 10 S	Gram-negative bacilli
Rapid 20E™	<i>Enterobacteriaceae</i>
API® 20 NE	Gram-negative non- <i>Enterobacteriaceae</i>
API® Staph	Staphylococci
API® 20 Strep	Streptococci
API® Candida	Yeasts
API® 20 C AUX	Yeasts
API® 20 A	Anaerobes
API® Coryne	Corynebacteria
API® Campy	<i>Campylobacter</i>
API® Listeria	<i>Listeria</i>
API® NH	<i>Neisseria, Haemophilus</i>
API® 50 CHE	<i>Enterobacteriaceae</i>
API® 50 CHL	Lactic bacteria
API® 50 CHB	<i>Bacillus</i>
ID 32 E	Gram-negative bacilli
Rapid ID 32 E	<i>Enterobacteriaceae</i>
ID 32 GN	Gram-negative bacilli
ID 32 STAPH	Staphylococci
Rapid ID 32 STREP	Streptococci
ID 32 C	Yeasts
Rapid ID 32 A	Anaerobes

- Species present in database
- Species identifiable with additional tests
- Changed taxon
- New species

	API® 20 E	API® 10 S	Rapid 20 E™	API® 20 NE			API® Staph	API® 20 Strep	API® Candida	API® 20 C AUX	API® 20 A	API® Coryne	API® Campy	API® Listeria	API® NH	API® 50 CHE	API® 50 CHL	API® 50 CHB		ID 32 E	Rapid ID 32 E	ID 32 GN	ID 32 Staph	Rapid ID 32 Strep	ID 32 C	Rapid ID 32 A			
<i>Abiotrophia defectiva</i>																													
<i>Achromobacter denitrificans</i>																													
<i>Achromobacter xylosoxidans</i>																													
<i>Acinetobacter/Moraxella spp</i>																													
<i>Acinetobacter/Pseudomonas spp</i>																													
<i>Acinetobacter baumannii</i>																													
<i>Acinetobacter baumannii/calcoaceticus</i>																													
<i>Acinetobacter calcoaceticus</i>																													
<i>Acinetobacter haemolyticus</i>																													
<i>Acinetobacter johnsonii</i>																													
<i>Acinetobacter junii</i>																													
<i>Acinetobacter junii/johnsonii</i>																													
<i>Acinetobacter lwoffii</i>																													
<i>Acinetobacter radioresistens</i>																													
<i>Acinetobacter spp</i>																													
<i>Actinobacillus pleuropneumoniae</i>																													
<i>Actinomyces israelii</i>																													
<i>Actinomyces meyeri</i>																													
<i>Actinomyces meyeri/odontolyticus</i>																													
<i>Actinomyces naeslundii</i>																													
<i>Actinomyces neuii ssp anitratus</i>																													
<i>Actinomyces neuii ssp neuii</i>																													
<i>Actinomyces odontolyticus</i>																													
<i>Actinomyces radingae</i>																													
<i>Actinomyces turicensis</i>																													
<i>Actinomyces viscosus</i>																													
<i>Aeromonas caviae</i>																													
<i>Aeromonas hydrophila/caviae</i>																													
<i>Aeromonas hydrophila</i>																													
<i>Aeromonas hydrophila/caviae/sobria</i>																													
<i>Aeromonas salmonicida ssp achromogenes</i>																													
<i>Aeromonas salmonicida masoucida/achromogenes</i>																													
<i>Aeromonas salmonicida ssp masoucida</i>																													
<i>Aeromonas salmonicida ssp salmonicida</i>																													
<i>Aeromonas sobria</i>																													
<i>Aerococcus urinae</i>																													
<i>Aerococcus viridans</i>																													
<i>Aggregatibacter actinomycetemcomitans</i>																													
<i>Aggregatibacter aphrophilus</i>																													
<i>Alcaligenes faecalis</i>																													
<i>Alcaligenes spp</i>																													
<i>Alloiooccus otitis</i>																													
<i>Anaerobiospirillum succiniciproducens</i>																													
<i>Anaerococcus prevotii</i>																													
<i>Aneurinibacillus aneurinilyticus</i>																													
<i>Arcanobacterium haemolyticum</i>																													
<i>Arcobacter cryaerophilus</i>																													
<i>Arthrobacter spp</i>																													
<i>Avibacterium paragallinarum</i>																													
<i>Bacillus amyloliquefaciens</i>																													
<i>Bacillus anthracis</i>																													

	API® 20 E	API® 10 S	Rapid 20 E™	API® 20 NE			API® Staph	API® 20 Strep	API® Candida	API® 20 C AUX	API® 20 A	API® Coryne	API® Campy	API® Listeria	API® NH	API® 50 CHE	API® 50 CHL	API® 50 CHB		ID 32 E	Rapid ID 32 E	ID 32 GN	ID 32 Staph	Rapid ID 32 Strep	ID 32 C	Rapid ID 32 A		
<i>Bacillus badius</i>																												
<i>Bacillus cereus</i>																												
<i>Bacillus circulans</i>																												
<i>Bacillus coagulans</i>																												
<i>Bacillus firmus</i>																												
<i>Bacillus lentus</i>																												
<i>Bacillus licheniformis</i>																												
<i>Bacillus megaterium</i>																												
<i>Bacillus mycoides</i>																												
<i>Bacillus pumilus</i>																												
<i>Bacillus smithii</i>																												
<i>Bacillus subtilis</i>																												
<i>Bacillus subtilis/amyloliquefaciens</i>																												
<i>Bacillus thuringiensis</i>																												
<i>Bacteroides caccae</i>																												
<i>Bacteroides eggerthii</i>																												
<i>Bacteroides fragilis</i>																												
<i>Bacteroides ovatus</i>																												
<i>Bacteroides ovatus/thetaiotaomicron</i>																												
<i>Bacteroides stercoris/eggerthii</i>																												
<i>Bacteroides stercoris</i>																												
<i>Bacteroides thetaiotaomicron</i>																												
<i>Bacteroides uniformis</i>																												
<i>Bacteroides vulgatus</i>																												
<i>Bacillus non reactive</i>																												
<i>Bergeyella zoohelcum</i>																												
<i>Bibersteinia trehalosi</i>																												
<i>Bifidobacterium adolescentis</i>																												
<i>Bifidobacterium bifidum</i>																												
<i>Bifidobacterium breve</i>																												
<i>Bifidobacterium dentium</i>																												
<i>Bifidobacterium longum</i>																												
<i>Bifidobacterium spp</i>																												
<i>Bordetella/Alcaligenes/Moraxella spp</i>																												
<i>Bordetella avium</i>																												
<i>Bordetella bronchiseptica</i>																												
<i>Bordetella spp</i>		</																										

- Species present in database
- Species identifiable with additional tests
- Changed taxon
- New species

	API® 20 E	API® 10 S	Rapid 20 E™	API® 20 NE		API® Staph	API® 20 Strep	API® Candida	API® 20 C AUX	API® 20 A	API® Coryne	API® Campy	API® Listeria	API® NH	API® 50 CHE	API® 50 CHL	API® 50 CHB		ID 32 E	Rapid ID 32 E	ID 32 GN	ID 32 Staph	Rapid ID 32 Strep	ID 32 C	Rapid ID 32 A
<i>Streptococcus agalactiae</i>							●																	●	
<i>Streptococcus alactolyticus</i>																								●	
<i>Streptococcus anginosus</i>							●																	●	
<i>Streptococcus bovis</i>																								●	
<i>Streptococcus canis</i>							●																	●	
<i>Streptococcus constellatus ssp constellatus</i>							●																	●	
<i>Streptococcus constellatus ssp pharyngis</i>							●																	●	
<i>Streptococcus constellatus</i>										●														●	
<i>Streptococcus downei</i>											●														●
<i>Streptococcus downei/sobrinus</i>																									●
<i>Streptococcus dysgalactiae ssp dysgalactiae</i>							●																	●	
<i>Streptococcus dysgalactiae ssp equisimilis</i>							●																	●	
<i>Streptococcus equi ssp equi</i>							●																	●	
<i>Streptococcus equi ssp zooepidemicus</i>							●																	●	
<i>Streptococcus equinus</i>							●																	●	
<i>Streptococcus gallolyticus ssp gallolyticus</i>							●																	●	
<i>Streptococcus gallolyticus ssp pasteurianus</i>							●																	●	
<i>Streptococcus gordonii</i>																									●
<i>Streptococcus group L</i>							●																		●
<i>Streptococcus infantarius ssp coli (Str.lutetiensis)</i>							●																	●	
<i>Streptococcus infantarius ssp infantarius</i>							●																	●	
<i>Streptococcus intermedius</i>							●				●													●	
<i>Streptococcus mitis</i>							●																●	●	
<i>Streptococcus mutans</i>							●																●	●	
<i>Streptococcus oralis</i>							●																●	●	
<i>Streptococcus parasanguinis</i>																								●	
<i>Streptococcus plurimalium</i>							●																	●	
<i>Streptococcus pneumoniae</i>							●																	●	
<i>Streptococcus porcinus</i>							●																	●	
<i>Streptococcus pyogenes</i>							●																	●	
<i>Streptococcus salivarius ssp salivarius</i>							●																	●	
<i>Streptococcus salivarius ssp thermophilus</i>							●								●									●	
<i>Streptococcus sanguinis</i>							●																	●	
<i>Streptococcus sobrinus</i>							●																	●	
<i>Streptococcus suis I</i>							●																	●	
<i>Streptococcus suis II</i>							●																	●	
<i>Streptococcus uberis</i>							●																	●	
<i>Streptococcus vestibularis</i>																								●	
<i>Tatumella ptyseos</i>																								●	
<i>Tetragenococcus halophilus</i>																●									●
<i>Trichosporon asahii</i>										●															●
<i>Trichosporon asteroides</i>																								●	
<i>Trichosporon inkin</i>										●														●	
<i>Trichosporon mucoides</i>										●														●	
<i>Trichosporon ovoides</i>																								●	
<i>Trichosporon spp</i>								●																	●
<i>Trueperella bernardiae</i>											●														●
<i>Trueperella pyogenes</i>											●														●
<i>Turicella otitidis</i>											●														●
<i>Vibrio alginolyticus/parahaemolyticus</i>															●										●
<i>Vibrio alginolyticus</i>	●		●	●											●						●	●			

	API® 20 E	API® 10 S	Rapid 20 E™	API® 20 NE		API® Staph	API® 20 Strep	API® Candida	API® 20 C AUX	API® 20 A	API® Coryne	API® Campy	API® Listeria	API® NH	API® 50 CHE	API® 50 CHL	API® 50 CHB		ID 32 E	Rapid ID 32 E	ID 32 GN	ID 32 Staph	Rapid ID 32 Strep	ID 32 C	Rapid ID 32 A
<i>Vibrio cholerae</i>	●		●	●											●						●	●			
<i>Vibrio fluvialis</i>	●		●	●											●						●	●			
<i>Vibrio metschnikovii</i>															●						●	●			
<i>Vibrio mimicus</i>	●			●											●						●	●			
<i>Vibrio parahaemolyticus</i>	●														●						●	●			
<i>Vibrio vulnificus/cholerae</i>			●												●						●	●			
<i>Vibrio vulnificus</i>	●			●											●						●	●			
<i>Veillonella parvula</i>										●															
<i>Veillonella spp</i>																									●
<i>Virgibacillus pantothenicus</i>																	●								
<i>Weeksella virosa</i>																								●	
<i>Weeksella virosa/Bergeyella zoohelcum</i>																								●	
<i>Weeksella virosa/Empedobacter brevis</i>				●																					
<i>Weissella confusa</i>																●									
<i>Weissella viridescens</i>																●									
<i>Yersinia aldovae</i>															●										
<i>Yersinia enterocolitica</i>	●	●	●												●						●	●			
<i>Yersinia frederiksenii/intermedia</i>	●														●						●	●			
<i>Yersinia frederiksenii</i>	●														●						●	●			
<i>Yersinia intermedia</i>	●														●						●	●			
<i>Yersinia kristensenii</i>	●														●						●	●			
<i>Yersinia pestis</i>	●														●						●	●			
<i>Yersinia pseudotuberculosis</i>	●	●	●												●						●	●			
<i>Yersinia ruckeri</i>	●														●						●	●			
<i>Zygosaccharomyces spp</i>																									●

YEASTS NOMENCLATURE

Name in the database	Imperfect State	Perfect State	Other names
<i>Candida albicans</i>	<i>Candida albicans</i>		
<i>Candida boidinii</i>	<i>Candida boidinii</i>		
<i>Candida catenulata</i>	<i>Candida catenulata</i>		
<i>Candida ciferrii</i>	<i>Candida ciferrii</i>	<i>Trichomonascus ciferrii</i>	<i>Stephanoascus ciferrii</i>
<i>Candida colliculosa</i>	<i>Candida colliculosa</i>	<i>Torulaspora delbrueckii</i>	
<i>Candida dattila</i>	<i>Candida dattila</i>	<i>Lachancea thermotolerans</i>	<i>Kluyveromyces thermotolerans</i>
<i>Candida dubliniensis</i>	<i>Candida dubliniensis</i>		
<i>Candida famata</i>	<i>Candida famata</i>	<i>Debaryomyces hansenii</i>	
<i>Candida glabrata</i>	<i>Candida glabrata</i>		
<i>Candida globosa</i>	<i>Candida globosa</i>	<i>Citeromyces matritensis</i>	
<i>Candida guilliermondii</i>	<i>Candida guilliermondii</i>	<i>Meyerozyma guilliermondii</i>	
<i>Candida hellenica</i>	<i>Candida hellenica</i>	<i>Zygoascus meyeræ</i>	
<i>Candida holmii</i>	<i>Candida holmii</i>	<i>Kazachstania exigua</i>	<i>Saccharomyces exiguus</i>
<i>Candida inconspicua</i>	<i>Candida inconspicua</i>		
<i>Candida intermedia</i>	<i>Candida intermedia</i>		
<i>Candida kefyr</i>	<i>Candida kefyr</i>	<i>Kluyveromyces marxianus</i>	
<i>Candida krusei</i>	<i>Candida krusei</i>	<i>Pichia kudriavzevii</i>	
<i>Candida lambica</i>	<i>Candida lambica</i>	<i>Pichia fermentans</i>	
<i>Candida lipolytica</i>	<i>Candida lipolytica</i>	<i>Yarrowia lipolytica</i>	
<i>Candida lusitanae</i>	<i>Candida lusitanae</i>	<i>Clavispora lusitanae</i>	
<i>Candida magnoliae</i>	<i>Candida magnoliae</i>		
<i>Candida melibiosica</i>	<i>Candida melibiosica</i>		
<i>Candida membranifaciens</i>	<i>Candida membranifaciens</i>		
<i>Candida norvegensis</i>	<i>Candida norvegensis</i>	<i>Pichia norvegensis</i>	
<i>Candida norvegica</i>	<i>Candida norvegica</i>		
<i>Candida parapsilosis</i>	<i>Candida parapsilosis</i>		
<i>Candida pelliculosa</i>	<i>Candida pelliculosa</i>	<i>Wickerhamomyces anomalus</i>	
<i>Candida pulcherrima</i>	<i>Candida pulcherrima</i>	<i>Metschnikowia pulcherrima</i>	
<i>Candida rugosa</i>	<i>Candida rugosa</i>		
<i>Candida sake</i>	<i>Candida sake</i>		

Name in the database	Imperfect State	Perfect State	Other names
<i>Candida silvicola</i>	<i>Candida silvicola</i>	<i>Nakazawaea holstii</i>	
<i>Candida spherica</i>	<i>Candida spherica</i>	<i>Kluyveromyces lactis</i> var <i>lactis</i>	
<i>Candida thermophila</i>	<i>Candida thermophila</i>	<i>Ogataea polymorpha</i>	<i>Hansenula polymorpha</i> , <i>Pichia angusta</i>
<i>Candida tropicalis</i>	<i>Candida tropicalis</i>		
<i>Candida utilis</i>	<i>Candida utilis</i>	<i>Lindnera jadinii</i>	
<i>Candida valida</i>	<i>Candida valida</i>	<i>Pichia membranifaciens</i>	
<i>Candida zeylanoides</i>	<i>Candida zeylanoides</i>		
<i>Cryptococcus albidus</i>	<i>Cryptococcus albidus</i>		
<i>Cryptococcus curvatus</i>	<i>Cryptococcus curvatus</i>		
<i>Cryptococcus humicola</i>	<i>Cryptococcus humicola</i>		
<i>Cryptococcus laurentii</i>	<i>Cryptococcus laurentii</i>		
<i>Cryptococcus neoformans</i>	<i>Cryptococcus neoformans</i>	<i>Filobasidiella neoformans</i>	
<i>Cryptococcus terreus</i>	<i>Cryptococcus terreus</i>		
<i>Cryptococcus uniguttulatus</i>	<i>Cryptococcus uniguttulatus</i>	<i>Filobasidium uniguttulatus</i>	
<i>Geotrichum candidum</i>	<i>Geotrichum candidum</i>	<i>Galactomyces candidus</i>	
<i>Geotrichum fermentans</i>	<i>Geotrichum fermentans</i>		
<i>Geotrichum klebahnii</i>	<i>Geotrichum klebahnii</i>		
<i>Kloeckera apiculata</i>	<i>Kloeckera apiculata</i>	<i>Hanseniaspora uvarum</i>	
<i>Kloeckera apis</i>	<i>Kloeckera apis</i>	<i>Hanseniaspora guilliermondii</i>	
<i>Kloeckera japonica</i>	<i>Kloeckera japonica</i>	<i>Hanseniaspora valbyensis</i>	
<i>Kodamaea ohmeri</i>	<i>Kodamaea ohmeri</i>		
<i>Lachancea kluyverii</i>		<i>Lachancea kluyverii</i>	<i>Saccharomyces kluyveri</i>
<i>Lindnera saturnus</i>		<i>Lindnera saturnus</i>	
<i>Milleromyza farinosa</i>	<i>Milleromyza farinosa</i>		
<i>Priceomyces carsonii</i>		<i>Priceomyces carsonii</i>	<i>Pichia carsonii</i>
<i>Prototheca wickerhamii</i>	<i>Prototheca wickerhamii</i>		
<i>Rhodotorula glutinis</i>	<i>Rhodotorula glutinis</i>		
<i>Rhodotorula minuta</i>	<i>Rhodotorula minuta</i>		
<i>Rhodotorula mucilaginosa</i>	<i>Rhodotorula mucilaginosa</i>		
<i>Saccharomyces cerevisiae</i>		<i>Saccharomyces cerevisiae</i>	
<i>Saprochaete capitata</i>	<i>Saprochaete capitata</i>	<i>Magnusiomyces capitatus</i>	
<i>Schwanniomyces etchellsii</i>		<i>Schwanniomyces etchellsii</i>	<i>Pichia etchellsii</i>
<i>Schwanniomyces polymorphus</i>		<i>Schwanniomyces polymorphus</i>	<i>Debaryomyces polymorphus</i>
<i>Sporobolomyces salmonicolor</i>	<i>Sporobolomyces salmonicolor</i>	<i>Sporidiobolus salmonicolor</i>	
<i>Trichosporon asahii</i>	<i>Trichosporon asahii</i>		
<i>Trichosporon asteroides</i>	<i>Trichosporon asteroides</i>		
<i>Trichosporon inkin</i>	<i>Trichosporon inkin</i>		
<i>Trichosporon mucoides</i>	<i>Trichosporon mucoides</i>		
<i>Trichosporon ovoides</i>	<i>Trichosporon ovoides</i>		

REAGENTS TO BE ORDERED

ADDITIONAL REAGENTS	API 20 E* (a)		Rapid ID 20E		API 10 S		API 20 NE		API Staph		API 20 Strep		API Coryne		API Listeria		API Candida		API 20 C AUX		API 20 A®		API Campy		API NH		API ZYM		API 50 CH		API 50 CH/E		API 20 E (a)		API 50 CH 50 CHL		
	20100	20700	10100	20050	20500	20600	20900	10300	10500	20210	20300	20800	10400	25200	50300-50430	20100	50300-50410																				
Zn (2x10g)	X			X																																	
TDA (2x1 ampule)	1/8		1/4																																		
VP 1 VP 2 (2x2 ampules)	1/8	1/8			1/8	1/8		1/20																													
NIT 1 NIT 2 (2x2 ampules)	1/8		1/4	1/8	1/8		1/20																														
ZYM A (2x1 ampule)					1/8	1/8	1/20						1/20																								
ZYM B (2x1 ampule)					1/8	1/8	1/20																														
NIN (2 ampules)							1/4																														
BCP (1 ampule)																1/2	1/4																				
EHR (1 ampule)																1/2																					
XYL (2 ampules)																1/8																					
James (2x1 ampule)	1/8	1/8	1/4	1/8																																	
FB (2x1 ampule)																																					
PYZ (2 ampules)																1/20	X																				
SUSPENSION MEDIUM																																					
NaCl 0.85% Medium (5 ml)																																					
NaCl 0.85% Medium (3 ml)																																					
NaCl 0.85% Medium (2 ml)			1/4		1/4																																
Suspension Medium (5 ml)	1/4			1/2																							1/4									1/10	
Suspension Medium (2 ml)																																				1/10	
ADDITIONAL PRODUCTS																																					
Oxidase Reagent	X	X	X	X																																	
API® OF Medium					X																																
API® M Medium					X																																
Paraffin oil	X	X	X	X	X	X	X	X		X		X	X	X													X									X	
Rack of 12 ampules	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Sterile PSipettes (5ml)	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Swabs																																					
Mc Farland standard		X		X	X	X				X	X	X	X	X																							

*figures = recommended quantity of additional reagent(s) for each strip kit ordered
(eg. TDA = 1/8 ↔ 1/8 of kit ref. 70402 (2 TDA ampoules) is required to use one kit of ref. 20100 (25 API 20 E strips))

X = additional product required
(a) : for the API 20 E strip, use the additional reagent kit ref. 20120 (7 ampules)

REAGENTS TO BE ORDERED

ADDITIONAL REAGENTS	ID 32 Staph		ID 32 C		ID 32 E		Rapid ID 32 E		Rapid ID 32 A		Rapid ID 32 Strep	
	32500	32200	32400	32700	32300	32600						
VP A VP B (2x2 ampules)	1/8											1/8
NIT 1 NIT 2	1/8						1/8					
NIN												1/8
James				1/8	1/8		1/8					
FB	1/8						1/8					1/8

SIMPLIFIED METHODOLOGIES

STRIPS	RÉF	MICRO ORGANISMS	SUSPENSION	MC FARLAND	TRANSFERT	MEDIUM	INCUBATION	ATMOSPHERE
API® 20 E	20100 20160	Enterobacteriaceae and other non-fastidious Gram negative bacilli	NaCl 0.85% Medium 5 ml ou Suspension Medium 5 ml	1 colony	NA	NA	37°C 18-24h/48h	Aerobic
RapiD 20 E	20701	Enterobacteriaceae	NaCl 0.85% Medium 2 ml	0.5 McF	NA	NA	37°C 4h	Aerobic
API® 10 S	10100	Enterobacteriaceae and other non-fastidious Gram negative bacilli	NaCl 0.85% Medium 5 ml ou Suspension Medium 5 ml	1 colony	NA	NA	37°C 18-24h	Aerobic
API® 20 NE	20050	Non-Enterobacteriaceae and non-fastidious Gram negative bacilli	NaCl 0.85% Medium 2 ml	0.5 McF	200 µL	API AUX	30°C 24-48h	Aerobic
API® Staph	20500	Genres Staphylococcus, Kocuria and Micrococcus	API Staph Medium	0.5 McF	NA	Medium	37°C 18-24h	Aerobic
API® 20 Strep	20600	Streptocoques	Suspension Medium 2 ml	4 McF	500 µl	NA	37°C 4h-24h	Aerobic
API® Coryne	20900	Bacteries coryneformes	Suspension Medium 3 ml	>6 McF	500 µl	Medium	37°C 24h	Aerobic
API® Listeria	10300	Listeria	Suspension Medium 2 ml	1 McF	NA	API GP	37°C 24h	Aerobic
API® Candida	10500	Yeasts	NaCl 0.85% Medium 2 ml	3 McF	NA	Medium	37°C 18-24h	Aerobic
API® 20 C AUX	20210	Yeasts	NaCl 0.85% Medium 2 ml Suspension Medium 2 ml	2 McF	100 µl	NA	30°C 48-72h	Aerobic
API® 20 A	20300	Anaerobes	API 20 A Medium	3 McF	NA	NA	37°C 24h	Anaerobic
API® Campy	20800	Campylobacter	NaCl 0.85 % Medium 3 ml	6 McF	150 µl	API C	37°C 24h	Aerobic/CO ₂
API® NH	10400	Neisseria – Haemophilus and Moraxella (<i>Branhamella</i>) catarrhalis	NaCl 0.85% Medium 2 ml	4 McF	NA	Medium	37°C 2h	Aerobic
API® ZYM	25200	Semi-quantitative enzyme activity tests	NaCl 0.85% Medium 2 ml ou Suspension Medium 2 ml	5 6 McF	NA	NA	37°C 4h	Aerobic
API® 50 CH API® 50 CHB/E	50300 50430	Bacillus	NaCl 0.85% Medium 5 ml + API 50 CHB/E Medium	2 McF 2 McF	NA	API AUX	30°C 24-48h 55°C 3/6-24h	Aerobic
API® 50 CH API® 50 CHB/E	50300 50430	Enterobacteriaceae	NaCl 0.85% Medium 5 ml + API 50 CHB/E Medium	0.5 McF 0.5 McF	NA	Medium	37°C 24-48h	Aerobic
API® 50 CH API® 50 CHL	50300 50410	Lactobacillus	API 50 CHL Medium	2 McF	NA	NA	30/37°C 48h	Aerobic

SIMPLIFIED METHODOLOGIES

STRIPS	RÉF	BACTERIA	SUSPENSION INOCULATOR 3 ML MANUAL 2 ML	MC FARLAND	TRANSFERT	MEDIUM	VOLUME DISPENSED	TEMP. INCUBATION	ATMOSPHERE
Rapid ID 32 E	32700	Enterobacteriaceae	NaCl 0.85% Medium	0.5	NA	NA	55 µl x 32	37°C 4h	Aerobic
Rapid ID 32 A	20701	Enterobacteriaceae	Suspension Medium	0.5 McF	NA	NA	55 µl x 32	37°C 4h	Aerobic
ID 32 E	10100	Streptococci	Suspension Medium	1 colonie	NA	NA	55 µl x 32	37°C 4h	Aerobic
ID 32 STAPH	20050	Enterobacteriaceae Non Enterobacteriaceae	NaCl 0.85% Medium	0.5 McF	NA	NA	55 µl x 32	37°C 24h	Aerobic + humidity
API® Staph	20500	Staphylococci	Suspension Medium	0.5 McF	NA	NA	55 µl x 32	37°C 24h	Aerobic + humidity
ID 32 C	20600	Yeasts	Suspension Medium	4 McF	250 µl	API C Medium	135 µl x 32	30°C 24/48h	Aerobic + humidity

BIBLIOGRAPHY

1 • *Elizabethkingia meningoseptica*

KIM (K.K.), KIM (M.K.), LIM (J.H.), PARK (H.Y.) and LEE (S.T.): Transfer of *Chryseobacterium meningosepticum* and *Chryseobacterium miricola* to *Elizabethkingia* gen. nov. as *Elizabethkingia meningoseptica* comb. nov. and *Elizabethkingia miricola* comb. nov. Int. J. Syst. Evol. Microbiol., 2005, 55, 1287-1293.

2 • *Kluyvera intermedia*

PAVAN (M.E.), FRANCO (R.J.), RODRIGUEZ (J.M.), GADALETA (P.), ABBOTT (S.L.), JANDA (J.M.) and ZORZÓPULOS (J.): Phylogenetic relationships of the genus *Kluyvera*: transfer of *Enterobacter intermedium* Izard et al. 1980 to the genus *Kluyvera* as *Kluyvera intermedia* comb. nov. and reclassification of *Kluyvera cochleae* as a later synonym of *K. intermedia*. Int. J. Syst. Evol. Microbiol., 2005, 55, 437-442.

3 • *Salmonella* spp

TINDALL (B.J.), GRIMONT (P.A.D.), GARRITY (G.M.) and EUZÉBY (J.P.): Nomenclature and taxonomy of the genus *Salmonella*. Int. J. Syst. Evol. Microbiol., 2005, 55, 521-524.

4 • *Bibersteinia trehalosi*

BLACKALL (P.J.), BOJENSEN (A.M.), CHRISTENSEN (H.) and BISGAARD (M.): Reclassification of [*Pasteurella*] *trehalosi* as *Bibersteinia trehalosi* gen. nov., comb. nov. Int. J. Syst. Evol. Microbiol., 2007, 57, 666-674.

5 • *Cronobacter* spp

IVERSEN (C.), MULLANE (N.), McCARDLELL (B.), TALL (B.D.), LEHNER (A.), FANNING (S.), STEPHAN (R.) and JOOSTEN (H.): *Cronobacter* gen. nov., a new genus to accommodate the biogroups of *Enterobacter sakazakii*, and proposal of *Cronobacter sakazakii* gen. nov., comb. nov., *Cronobacter malonaticus* sp. nov., *Cronobacter turicensis* sp. nov., *Cronobacter muytjensii* sp. nov., *Cronobacter dublinensis* sp. nov., *Cronobacter genomospecies 1*, and of three subspecies, *Cronobacter dublinensis* subsp. *dublinensis* subsp. nov., *Cronobacter dublinensis* subsp. *lausannensis* subsp. nov. and *Cronobacter dublinensis* subsp. *lactaridi* subsp. nov. Int. J. Syst. Evol. Microbiol., 2008, 58, 1442-1447.

6 • *Streptococcus salivarius* ssp *salivarius*

HOWEY (R.T.), LOCK (C.M.) and MOORE (L.V.H.): Subspecies names automatically created by Rule 46. Int. J. Syst. Bacteriol., 1990, 40, 317-319.

7 • *Streptococcus salivarius* ssp *thermophilus*

VALIDATION LIST N° 54. Int. J. Syst. Bacteriol., 1995, 45, 619-620. [SCHLEIFER (K.H.), EHRMANN (M.), KRUSCH (U.) and NEVE (H.): Revival of the species *Streptococcus thermophilus* (ex *Orla-Jensen*, 1919) nom. rev. Syst. Appl. Microbiol., 1991, 14, 386-388.] Note: The Validation List N° 54 [1] encompasses the following citation: «Name: *Streptococcus thermophilus* (basonym: *Streptococcus salivarius* subsp. *thermophilus*) - Proposed as: Revived name - Authors: Schleifer et al. 1991 [2] - Nomenclatural type: ATCC 19258 (= NCTD 573)». However, *Streptococcus thermophilus* cannot be a revived name because it appeared in the Approved Lists (1980). Obviously, *Streptococcus thermophilus* nom. rev. (sic) is erroneously cited in the Validation List N° 54.

8 • *Streptococcus bovis*

Streptococcus gallolyticus ssp *gallolyticus* *Streptococcus gallolyticus* ssp *pasteurianus* *Streptococcus infantarius* ssp *coli* *Streptococcus infantarius* ssp *infantarius* SCHLEGEL (L.), GRIMONT (F.), COLLINS (M.D.), RÉGNAULT (B.), GRIMONT (P.A.D.) and BOUVET (A.): *Streptococcus infantarius* sp. nov., *Streptococcus infantarius* subsp. *infantarius* subsp. nov. and *Streptococcus infantarius* subsp. *coli* subsp. nov., isolated from humans and

food. Int. J. Syst. Evol. Microbiol., 2000, 50, 1425-1434. ASSOCIATE EDITOR, IJSB: Notification that new names and new combinations have appeared in volume 50, part 4, of the IJSEM. Int. J. Syst. Evol. Microbiol., 2000, 50, 1701-1702.

POYART (C.), QUESNE (G.) and TRIEU-CUOT (P.): Taxonomic dissection of the *Streptococcus bovis* group by analysis of manganese-dependent superoxide dismutase gene (*sodA*) sequences: reclassification of *Streptococcus infantarius* subsp. *coli* as *Streptococcus lutetiensis* sp. nov. and of *Streptococcus bovis* biotype II.2 as *Streptococcus pasteurianus* sp. nov. Int. J. Syst. Evol. Microbiol., 2002, 52, 1247-1255.

SCHLEGEL (L.), GRIMONT (F.), AGERON (E.), GRIMONT (P.A.D.) and BOUVET (A.): Reappraisal of the taxonomy of the *Streptococcus bovis*/*Streptococcus equinus* complex and related species: description of *Streptococcus gallolyticus* subsp. *gallolyticus* subsp. nov., *S. gallolyticus* subsp. *macedonicus* subsp. nov. and *S. gallolyticus* subsp. *pasteurianus* subsp. nov. Int. J. Syst. Evol. Microbiol., 2003, 53, 631-645.

9 • **Freydiere AM.**: Yeast identification in the clinical microbiology laboratory: phenotypical methods. Medical Mycology, 2001, Vol39, (1), p9-33

10 • **Budak A.**: Epidemiology of *Candida* infection. II. Application of biochemical methods for typing of *Candida albicans* strains.

Arch Immunol Ther Exp (Warsz).1990, Vol 38, (5-6), p369-77

11 • **C.P.Kurtzman**: The Yeast, a taxonomic study. 5ème Edition

12 • **Saprochaete capitata**: Book : The Yeasts, Fifth Edition: A Taxonomic Study by C.P. Kurtzman, J.W. Fell

13 • *Parabacteroides distasonis*

SAKAMOTO (M.): Reclassification of *Bacteroides distasonis*, *Bacteroides goldsteinii* and *Bacteroides merdae* as *Parabacteroides distasonis* gen. nov., comb. Nov., *Parabacteroides goldsteinii* comb. Nov. and *Parabacteroides merdae* comb; nov. - International Journal of Systematic and Evolutionary Microbiology, 2006,56,1599-1605

14 • *Campylobacter ureolyticus*

VANDAMME (P): Reclassification of *Bacteroides ureolyticus* as *Campylobacter ureolyticus* comb. nov., and emended description of the genus *Campylobacter* - International Journal of Systematic and Evolutionary Microbiology, 2010,60,2016-2022

15 • *Parvimonas micra*

TINDALL (B.J.) and EUZÉBY (J.P.): Proposal of *Parvimonas* gen. nov. and *Quatronicoccus* gen. nov. as replacements for the illegitimate, prokaryotic, generic names *Micromonas* Murdoch and Shah 2000 and *Quadricoccus* Maszenan et al. 2002, respectively. Int. J. Syst. Evol. Microbiol., 2006, 56, 2711-2713.

16 • *Propionibacterium propionicum*

MOORE (W.E.C.) and MOORE (L.V.H.): Index of the bacterial and yeast nomenclatural changes published in the International Journal of Systematic Bacteriology since the 1980 Approved Lists of bacterial names (1 January 1980 to 1 January 1992). American Society for Microbiology, Washington, D.C., 1992.

17 • **Capnocytophaga species: FRANDSEN (EVG):** Diversity of *Capnocytophaga* species in children and description of *Capnocytophaga leadbetteri* sp. Nov. and *Capnocytophaga* genospecies AHN8471 - International Journal of Systematic and Evolutionary Microbiology, 2008,58,324-336. GILLIGAN (PH): *Capnocytophaga ochracea* Septicemia - Journal of Clinical Microbiology, 1981,13,4,643-645. WINN (RE): Septic arthritis involving *Capnocytophaga*

ochracea - Journal of Clinical Microbiology, 1984,19,4,538-540. HAWKEY (PM): *Capnocytophaga ochracea* infection: two cases and a review of the published work - Journal Clin Pathol, 1984,37,1066-1070. TANNER (ACR): API® ZYM and API® an-ident reactions of fastidious oral gram-negative species - Journal of Clinical Microbiology, 1985,22,3,333-335. RUMMENS (JL): Isolation of *Capnocytophaga* species with a new selective medium- Journal of Clinical Microbiology, 1985,22,3,375-378. BUU-HOI (AY): Endocarditis caused by *Capnocytophaga ochracea* - Journal of Clinical Microbiology, 1988,26,5,1061-1062. YAMAMOTO (T): *Capnocytophaga haemolytic* asp. Nov. and *Capnocytophaga granulosa* sp. Nov., from human dental plaque - International Journal of Systematic Bacteriology, 1994,44,2,324-329.

18 • **YASSIN (A.F.), HUPFER (H.), SIERING (C.) and SCHUMANN (P.):** Comparative chemotaxonomic and phylogenetic studies on the genus *Arcanobacterium* Collins et al. 1982 emend. Lehnen et al. 2006: proposal for *Trueperella* gen. nov. and emended description of the genus *Arcanobacterium*. Int. J. Syst. Evol. Microbiol., 2011, 61, 1265-1274.

19 • *Campylobacter sputorum*

ON (S.L.W.), ATABAY (H.I.), CORRY (J.E.L.), HARRINGTON (C.S.) and VANDAMME (P.): Emended description of *Campylobacter sputorum* and revision of its infrasubspecific (biovar) divisions, including *C. sputorum* biovar *paraureolyticus*, a urease-producing variant from cattle and humans. Int. J. Syst. Bacteriol., 1998, 48, 195-206. VANDAMME (P.) and ON (S.L.W.): Recommendations of the Subcommittee on the taxonomy of *Campylobacter* and related bacteria. Int. J. Syst. Evol. Microbiol. 2001, 51, 719-721.

20 • *Avibacterium paragallinarum*

BLACKALL (P.J.): Reclassification of *Pasteurella gallinarum*, [*Haemophilus*] *paragallinarum*, *Pasteurella avium* and *Pasteurella volantium* as *Avibacterium gallinarum* gen. nov., comb. nov., *Avibacterium paragallinarum* comb. nov., *Avibacterium avium* comb. nov. and *Avibacterium volantium* comb. nov. - International Journal of Systematic and Evolutionary Microbiology, 2005,55,353-62

21 • *Aggregatibacter actinomycetemcomitans*, *Aggregatibacter aphrophilus*

NORSKOV-LAURITSEN (N): Reclassification of *Actinobacillus actinomycetemcomitans*, *Haemophilus aphrophilus*, *Haemophilus paraprofundus* and *Haemophilus segnis* as *Aggregatibacter actinomycetemcomitans* gen. nov., comb. nov., *Aggregatibacter aphrophilus* comb. nov. and *Aggregatibacter segnis* comb. nov., and emended description of *Aggregatibacter aphrophilus* to include V factor-dependent and V factor-independent isolates - International Journal of Systematic and Evolutionary Microbiology, 2006,56,2135-46.

22 • *Lysinibacillus fusiformis* and *Lysinibacillus sphaericus*

AHMED (I), YOKOTA (A): Proposal of *Lysinibacillus boronitolerans* gen. nov. sp. nov., and transfer of *Bacillus fusiformis* to *Lysinibacillus fusiformis* comb. nov. and *Bacillus sphaericus* to *Lysinibacillus sphaericus* comb. nov. Int J Syst Evol Microbiol. 2007 May;57(Pt 5):1117-25.

23 • *Streptococcus salivarius*

FARROW (JA), COLLINS (MD): DNA base composition, DNA-DNA homology and long-chain fatty acid studies on *Streptococcus thermophilus* and *Streptococcus salivarius*. J Gen Microbiol. 1984 Feb;130(2):357-62.

24 • *Lactobacillus curvatus*

TORRIANI (S), VAN REENEN (CA), KLEIN (G): *Lactobacillus curvatus* subsp. *curvatus* subsp. nov. and *Lactobacillus curvatus* subsp. *melibiosus* subsp. nov. and *Lactobacillus sake* subsp. *sake* subsp. nov. and *Lactobacillus sake* subsp. *carnosus* subsp. nov., New Subspecies of *Lactobacillus curvatus* Abo-Elnaga and Kandler 1965 and *Lactobacillus sake* Katagiri, Kitahara, and Fukami 1934 (Klein et al. 1996, Emended Descriptions), Respectively. Int. J. Syst. Evol. Microbiol., 1996, 46, 4, 1158-1163

25 • *Pectobacterium* species :

GARDAN (L): Elevation of three subspecies of *Pectobacterium carotovorum* to species level: *Pectobacterium atrosepticum* sp. nov., *Pectobacterium betavascularum* sp. nov. and *Pectobacterium wasabia* sp. nov. - International Journal of Systematic and Evolutionary Microbiology, 2003, 53, p381-39.

26 • *SPROER (C) :*

The phylogenetic position of *Serratia*, *Buttiauxella* and some other genera of the family Entero-

bacteriaceae - International Journal of Systematic Bacteriology, 1999, 49, p1433-1438 HAUBEN (L) : Phylogenetic position of Phytopathogens within the Enterobacteriaceae System. Appl. Microbiol., 1998, 21, p384-397 MERGAERT (J) : Reclassification of non-pigmented *Erwinia herbicola* strains from trees as *Erwinia billingiae* sp. nov. - International Journal of Systematic Bacteriology, 1999, 49, p377- 383

27 • *Staphylococcus pseudintermedius* :

VAN HOOVELS (L): The first case of *Staphylococcus pseudintermedius* in humans isolated from an ICD lead - ESCMID, Nice, 2006, P1698 DEVRIESE (LA): *Staphylococcus pseudintermedius* sp. Nov., a coagulase-positive species from animals - IJSEM, 2005,55,1569-1573. SASAKI (T): Reclassification of phenotypically identified *Staphylococcus intermedius* strains - JCM, 2007,45,2770-2778. HAJEK (V): *Staphylococcus intermedius*, a new species isolated from animals - IJBS, 1976,26,401-408

28 • *Prototheca wickerhamii* :

PADHYE (A): Rapid identification of *Prototheca* species by API® 20C - Journal of Clinical Microbiology, 1979,10,4,579-582. McMULLAN (B): *Prototheca wickerhamii* mimicking yeast: a cautionary tale - Journal of Clinical Microbiology, 2011,49,8,3078-3081.

29 • *Parabacteroides distasonis* et *Parabacteroides merdae*

SAKAMOTO (M): Reclassification of *Bacteroides distasonis*, *Bacteroides goldsteinii* and *Bacteroides merdae* as *Parabacteroides distasonis* gen. nov., comb. Nov., *Parabacteroides goldsteinii* comb. Nov. and *Parabacteroides merdae* comb; nov. - International Journal of Systematic and Evolutionary Microbiology, 2006,56,1599-1605

30 • *Propionibacterium propionicum*

MOORE (W.E.C.) and MOORE (L.V.H.): Index of the bacterial and yeast nomenclatural changes published in the International Journal of Systematic Bacteriology since the 1980 Approved Lists of bacterial names (1 January 1980 to 1 January 1992). American Society for Microbiology, Washington, D.C., 1992

31 • *Pseudoflavonifractor capillosus*

CARLIER (J.P.), BEDORA-FAURE (M.), K'OUAS (G.), ALAUZET (C.) and MORY (F.): Proposal to unify *Clostridium orbiscindens* Winter et al. 1991 and *Eubacterium plautii* (Seguin 1928) Hofstad and Aasjord 1982, with description of *Flavonifractor plautii* gen. nov., comb. nov., and reassignment of *Bacteroides capillosus* to *Pseudoflavonifractor capillosus* gen. nov., comb. nov. Int. J. Syst. Evol. Microbiol., 2010, 60, 585-590.

32 • **LENNETTE E.H., BALOWS A., HAUSLER (JR) W. J., SHADOMY H.J.** Manual of Clinical Microbiology Fifth Edition (1991) Amer. Soc. Microbiol., Washington, D. C. 2. BARNETT J.H., PAYNE R. W., YARROW D. Yeasts: Characteristics and Identification (First edition 1983, Second 1990, Third 2000) Cambridge University Press - London

33 • **BARRETT T. J., PATTON C.M., MORRIS G.K**

Differentiation of *Campylobacter* species Using Phenotypic Characterization (1988) Laboratory Medicine, 19, 96-102

34 • **BOUVET P.J.M., GRIMONT P.A.D.** Taxonomy of the Genus *Acinetobacter* with the Recognition of *Acinetobacter baumannii* sp. Nov., *Acinetobacter haemophilus* sp. Nov., *Acinetobacter johnsonii* sp. nov., and *Acinetobacter junii* sp. Nov. And Amended Descriptions of *Acinetobacter calcoaceticus* and *Acinetobacter lwoffii* (1986) Int. J. Syst. Bacteriol., 36, 228-240

35 • **FARMER J.J. III, DAVIS B.R., HICKMAN-BRENNER F.W., McWHORTER A., HUNTLEY-CARTER G.P., ASBURY M.A., WATHEN-GRADY H.G., ELIAS C., FANNING G.R., STEIGERWALT A.G., O'HARA C.M., MORRIS G.K., SMITH P.B., DON BRENNER J.** Biochemical Identification of New Species and Biogroups of Enterobacteriaceae Isolated from Clinical Specimens (1985) J. Clin. Microbiol., 21, 46-76

36 • **GILARDI G.L.**

Identification of Glucose Nonfermenting Gram Negative Rods (1989) - Document Interne 7.HOLDEMAN L.V., CATO E.P., MODRE W.E.C. Anaerobe Laboratory Manual 4th Edition (1977) Virginia Polytechnic Institute and State University - Virginia - USA

37 • HOLLIS D.G., WEAVER R.E.

Gram Positive organisms: a guide to identification

Atlanta Special Bacteriology Section(1984) Center for Disease Control

38 • HOLMES B., PINNING C.A., DAWSON C.A.

A Probability Matrix for the Identification of Gram-negative, Aerobic, Non-fermentative Bacteria that Grow on Nutrient agar (1986) J. Gen. Microbiol., 132, 1827-1842

39 • HOLMES B., PINNING C.A., DAWSON C.A.

A Revised Probability Matrix for the Identification of Gram-negative, Aerobic Rod-shaped, Fermentative Bacterial (1986) J. Gen. Microbiol., 132, 3113-3135

40 • KILIAN M., MIKKELSEN L., HENRICHSEN J.

Taxonomic Study of Viridans Streptococci: Description of *Streptococcus gordonii* sp. Nov. and Emended Descriptions of *Streptococcus sanguis* (White and Niven (1946)), *Streptococcus oralis* (Bridge and Sneath 1982), and *Streptococcus mitis* (Andrews and Horder 1906).(1986) Int. J. Bacteriol., 39, 471-484

41 • KRIEG N.R., HOLT J.G.

Bergey's Manual of Systematic Bacteriology Ninth Edition - Volume 1 (1984) Williams & Wilkins Co., Baltimore, Md

42 • MEGRAUD F., CHEVRIER D., DESPLACES N., SEDALLIAN A., GUESDON J.L.

Urease-Positive Thermophilic Campylobacter Isolated from an Appendix and from Human Feces. (1988) J. Clin. Microbiol., 26, 1050-1051

43 • MOORE L.V.H., MOORE W.E.C.

Anaerobe Lab Manual Update (1991) Virginia Polytechnic Institute and State University. Blacksburg, Virginia, USA

44 • PELOUX Y,CANIAUX I.

Identification des corynebactéries et germes apparentés. Apport d'une nouvelle galerie d'identification : API® Coryne (1990) Feuillet de Biologie, XXI, 27-34

45 • RICHARD C.

Nouvelles Espèces de Enterobacteriaceae (1984) Bull. Institut Pasteur, 82, 255-277

46 • SCHLEIFER K.H., KRAUS J., DVORAK C., KILPER-BAELTZ R., COLLINS M.D.

Transfer of *Streptococcus lactis* and Related Streptococci to the Genus *Lactococcus* gen. Nov. (1985) System. Appl. Microbiol., 6, 183-195

47 • SKINNER F.A., QUESNEL L.B.

Streptococci (1978) Soc. Appl. Bacteriol. Symposium Série N° 7 - Academic Press

48 • SNEATH P.H.A., MAIR N.S., SHARPE M.E., HOLT J.G.

Bergey's Manual of Systematic Bacteriology Ninth Edition - Volume 2 (1986) Williams & Wilkins Co., Baltimore, Md

49 • STEELE T. W., OWEN R.J.

Campylobacter jejuni subsp. doylei subsp. nov., a Subspecies of Nitrate-Negative Campylobacters Isolated from Human Clinical Specimens (1988) Int. J. Syst. Bacteriol., 38, 316-318

50 • KREGER VAN RIJ N. J. W.

«The Yeasts - A Taxonomic Study.» Livre Ed. ELSEVIER (1984)

51 • COYLE M.B., NOWOWIEJSKI D.J., RUSSEL J.Q., GROMAN N.B.

«Laboratory review of reference strains of *Corynebacterium diphtheriae* indicates mistyped intermedius strains.» (1993) J. Clin. Microbiol. 31, 11, 3060-3062.

52 • FUNKE G, LAWSON P.A, COLLINS M.D.

«Heterogeneity within human-derived centers for disease control and prevention (CDC) Coryneform group ANF- 1 - like bacteria and description of *Corynebacterium auris* sp. nov.» (1981) Int. J. Syst. Bact. 45, 4, 735-739.

53 • RIEGEL P., RUIMY R., DE BRIEL D., PREVOST G., JEHL F., CHRISTEN R., MONTEIL H.

«*Corynebacterium seminale* sp. nov., a new species of Coryneform organism associated with genital infection in male patients.» (1995) J. Clin. Microbiol. 33, 2244-2249.

54 • RIEGEL P., DE BRIEL D., PREVOST G., JEHL F., MONTEIL H., MINCK R. «Taxonomic Study of *Corynebacterium* group ANF-1 strains : Proposal of *Corynebacterium afermentans* sp. nov. containing the Subspecies *Corynebacterium afermentans* subsp. *afermentans* subsp. nov. and *Corynebacterium afermentans* subsp. nov. *lipophilum* subsp. nov.» (1993) Int. J. Syst. Bact. Vol.43, N°2, 287-292.

55 • RIEGEL P., DE BRIEL D., PREVOST G., JEHL F., MONTIEL H.

«Proposal of *Corynebacterium propinquum* sp. nov. for *Corynebacterium* group ANF-3 strains.» (1993) FEMS Microbiology Letters Vol. 113

56 • MURRAY P.R., BARON E.J., PFALLER M.A., TENOVER F.C., YOLKEN R.H. «Manual of clinical microbiology». (Sixth edition 1995 , 7th 1999 , 8th 2003 Livre ED. ASM PRESS.)

57 • RENAUD F.N.R.,DUTAUR M., DAOUD S., AUBEL D.,RIEGEL P., MONGET D., FRENEY J., Differentiation Of *Corynebacterium amycolatum*, *C. minutissimum* , and *C. striatum* by carbon assimilation tests. (1998) J.Clin.Microbiol., 36, 3698-370<Bibliographie>2.

58 • FUNKE G., VON GRAEVENITZ A., CLARRIDGE J.E., BERNARD K.A. «Clinical microbiology of Coryneform bacteria». (1997) Clinical microbiology reviews Vol10, 1, 125-159.

59 • DAOUD S.E

«Apport des tests auxanotrophiques dans l'identification de quatre espèces de Corynebactéries.» (1996) Diplôme d'études approfondies de Génie Biologique et Médical.

60 • R. A. WHILEY, H. Y. FRASER, C. W. I. DOUGLAS, J. M. HARDIE, A. M. WILLIAMS, and M. D. COLLINS

Streptococcus parasanguis sp. nov., an atypical viridans *Streptococcus* from human clinical specimens. FEMS Microbiology Letters 68 (1990) p. 115-121

61 • SHINJO (T.), FUJISAWA (T.) and MITSUOKA (T.):

Proposal of two subspecies of *Fusobacterium necrophorum* (Flügge) Moore and Holdeman: *Fusobacterium necrophorum* subsp. *necrophorum* subsp. nov., nom. rev. (ex Flügge 1886), and *Fusobacterium necrophorum* subsp. *funduliforme* subsp. nov., nom. rev. (ex Hallé 1898). Int. J. Syst. Bacteriol., 1991, 41, 395-397

62 • SCHLEGEL (L.), GRIMONT (F.), AGERON (E.), GRIMONT (P.A.D.) and BOUVET (A.):

Reappraisal of the taxonomy of the *Streptococcus bovis*/*Streptococcus equinus* complex and related species: description of *Streptococcus gallolyticus* subsp. *gallolyticus* subsp. nov., *S. gallolyticus* subsp. *macedonicus* subsp. nov. and *S. gallolyticus* subsp. *pasteurianus* subsp. nov. Int. J. Syst. Evol. Microbiol., 2003, 53, 631-645.

63 • POYART (C.), QUESNE (G.) and TRIEU-CUOT (P.): Taxonomic dissection of the *Streptococcus bovis* group by analysis of manganese-dependent superoxide dismutase gene (*sodA*) sequences: reclassification of *Streptococcus infantarius* subsp. *coli* as *Streptococcus lutetiensis* sp. nov. and of *Streptococcus bovis* biotype II.2 as *Streptococcus pasteurianus* sp. nov. Int. J. Syst. Evol. Microbiol., 2002, 52, 1247-1255.

64 • AGUIRRE (M.) and COLLINS (M.D.): Phylogenetic analysis of some *Aerococcus*-like organisms for urinary tract infections: description of *Aerococcus urinae* sp. nov. J. Gen. Microbiol., 1992, 138, 401-405.

65 • COLLINS (M.D.), AGUIRRE (M.), FACKLAM (R.R.), SHALLCROSS (J.) and WILLIAMS (A.M.): *Globicatella sanguis* gen. nov., sp. nov., a new gram-positive catalase-negative bacterium from human sources. J. Appl. Bacteriol. 1992, 73, 433-437.]

66 • G GAUTHIER, B LAFAY, R RUIMY, V BREITTMAYER, JL NICOLAS, M GAUTHIER, and R CHRISTEN Small-subunit rRNA sequences and whole DNA relatedness concur for the reassignment of *Pasteurella piscicida* (Snieszko et al.) Janssen and Surgalla to the genus *Photobacterium* as *Photobacterium damsela* subsp. *piscicida* comb. Nov Int. J. Syst. Bacteriol., Jan 1995; 45: 139 - 144.