

INTRODUCTION AND PURPOSE

In diagnosis of hepatitis C infection, updated assays for anti-HCV detection perform very high on positive samples from chronic carriers infected with any viral genotype. Differences in performance can be recorded on samples from patients undergoing the early moments of the anti-HCV response after an acute primary infection or having very low anti-HCV levels (i.e. patients with anti-HCV in serum after a past resolved HCV infection). In both cases, differences between tests originate mostly from detection ability for antibodies to NS3 and at a second level to core and NS4.


The VIDAS® automated system (bioMérieux) is suitable for routine, emergency or complementary testing, notably in infectious diseases diagnosis.



To complete its existing HIV-HAV-HBV panel, VIDAS® Anti-HCV has been developed. Following the assessment of its specificity performance, we performed a sensitivity evaluation of this prototype by comparison with CE-marked tests, notably on well-characterized positive clinical samples.

METHODS

VIDAS® Anti-HCV is a third-generation test using antigens corresponding to HCV core, NS3 and NS4 for the qualitative detection of anti-HCV IgG antibodies. The CE-marked tests used in this study are Architect (Abbott) and Vitros (Ortho-Clinical diagnostics). Two immunoblots, RIBA (Ortho-Clinical diagnostics) and LIA (Innogenetics), were used to characterize antibody nature of samples.



ASSAYS	Format	Detection	1-191	1027-1657	1658-1711	1712-1972	1973-2420	2421-3011
VIDAS (ELFA)	Indirect sandwich	IgG	peptide (44aa)	recombinant protein (1192-1457)	peptide (43aa)	peptide (29aa)		
ARCHITECT (CLIA)	Indirect sandwich	IgG+IgM	E. coli 1-150	HC43 1192-1457	Yeast rec C100-3 1569-1931			
VITROS (CLEIA)	Indirect sandwich	IgG	peptide (C22) 10-53	Yeast rec C200 1192-1931		Yeast rec 2054-2995		
RIBA HCV (Immuno BLOT)	Indirect sandwich	IgG	peptide (C22) 10-53	c33c 1192-1457	peptide (C100p) 1694-1735	peptide (S-1-1p) 1920-1935		Yeast rec 2054-2995
LIA HCV (Immuno BLOT)	Indirect sandwich	IgG	peptides (C1) 1-32 (C2) 31-74	recombinant 1188-1465	peptide 1696-1739	peptide 1916-1944		peptide 2263-2318

FIGURE 1: Formats and representative HCV proteins used in the different Anti-HCV assays

The specificity study was performed on:

> 4766 French blood donor samples (drawn from different plants in France), HCV negatives after being tested with molecular biology assay. The false positive sera were characterized, notably those linked to the NS3 recombinant protein used in VIDAS® Anti-HCV by a 2D-PAGE analysis: on Figure 2, see the 2 different profiles of interferences: detection of molecules fixed directly on the NS3 sequence (A) or on co-purified proteins with NS3 (B).




FIGURE 2: 2D-PAGE analysis

The sensitivity study was performed on a collection of samples:


> 402 positive samples collected from patients at different stages of HCV disease and covering the 6 genotypes (1 to 6) (Biomnis Laboratory).
> 97 anti-HCV well-characterized sera, notably on their antibody nature (anti-core, anti-NS3) from chronic or resolved past infection patients.

All of them have been tested in comparison with Architect® or VITROS.

RESULTS

1. SPECIFICITY:

Performed on 4766 negative samples, this analysis demonstrated a high specificity performance of the VIDAS® Anti-HCV assay (>99.7%) with an excellent discrimination: >98 % of negatives have index under 0.5



Analysis of false positive shows that most of interferences are linked to NS3 raw material (>84.5%) which can be due to interferences recognizing directly either NS3 seq or contaminant co-purified proteins.

Number of false positives on 4766 donors	PTB3	NAME	Test Values				interference identified	
			PTB3	PTB4	PTB5	PTB6		
13 (100% with S/CO<5)	PTB3	EFS D676	1,08	0,65	0,82	0,37	NS3 (ND)	
	PTB3	EFS B0105	1,09	0,57	empty	empty	NS3 (ND)	
	PTB3	EFS B0267	1,15	0,66	empty	empty	NS3 (ND)	
	PTB3	S16-M005	1,18	0,73	0,56	0,32	NS3 E coli	
	PTB3	S11-M018	1,43	0,59	1,36	0,87	NS3 seq	
	PTB3	EFS B0114	1,88	1,88	1,04	1,11	0,71	NS3 (ND)
	PTB3	EFS B0274	2,40	2,4	2,36	1,59	2,47	NS3 seq
	PTB3	S11-F018	2,74	2,74	1,65	0,98	0,39	NS3 E coli
	PTB5	ECH1017	1,01	1,37	1,14	1,01	1,18	NS3 (ND)
	PTB5	EFS0311	1,17	1,34	1,44	1,17	1,21	NS3 seq
	PTB5	EFS0654	1,33	1,97	1,95	1,33	0,95	Format
	PTB5	EFS0291	2,48	4,62	4,66	2,48	4,76	NS3 seq
	PTB5	EFS0308	3,13	3,13	2,45	3,13	4,67	N27T

11/13 are linked to NS3

2. SENSITIVITY:

VIDAS® Anti-HCV sensitivity demonstrated an efficient antibody detection for the six HCV genotypes whatever the antibody titer and response,

VIDAS Anti-HCV	POSITIVE SERA			
	STUDY	TOTAL	LOW*	HIGH
	Positive samples	402 (100%)	52 (12,94%)	350 (87,06%)
% sensitivity	100%	100%	100%	
% samples with s/CO>5	92%	35%	100%	

* Low positive as CDC guidelines, i.e. <5 for Architect

with notably a performing anti-NS3 detection evidenced on the panel of the 97 anti-HCV well-characterized sera

SAMPLES	PCR	Immunoassays		BLOTS	
		VIDAS	VITROS	RIBA	LIA
anti-NS3 alone	-	25/30	17/29	27/30	24/30
anti-CORE alone	-	23/29	24/29	23/28	23/28
Weak positive (only two antibodies)	-	14*/16	16/16	5/16	16/16
Weak positive (only two antibodies)	+	3/3	3/3	2/3	2/2
Chronic infections	+	19/19	18/18	ND	ND
Total	75 RNA - 22 RNA +	84/97 86,6%	78/95 82,1%	57/77 74,0%	65/76 85,5%

* 2 samples have index within the grey zone (0,8-0,99)

CONCLUSION

Evaluation of VIDAS® Anti-HCV sensitivity and specificity showed that this new prototype VIDAS assay is as performant as other CE-marked tests and suitable for the qualitative detection of antibodies anti-HCV in human samples.